



**RAPID MR International, LLC**  
 Columbus, Ohio 43204, USA  
 1500 Lake Shore Drive, Suite 310  
 www.rapidmri.com

## Determining your Magnet's *Polarity*

Using a compass, please identify either 'QSN' or 'QNS' and communicate the result.

Thank you very much!

*Service end: coil/gradient cables and coolant hoses lead out here*  
*Patient end: examination subject is put into magnet from here*

*Hint: compare earth magnetic field  
 (geographic north pole (compass) = magnetic south pole)*

